

### **QuiKResponse - 1 liter Oil & Grease, waste water (WasteWatR TM)**

Art. ID                      ERA-518QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

#### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One liter whole-volume glass bottle with a 33-430 thread is ready to analyze. For SPE compatible QC (33-400 thread), use Cat #818 1 liter Boston Round Oil & Grease. Please specify if for SPE compatible PT.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Oil and Grease		20 - 100	mg/L		