

## **QuiKResponse - HEM / SGT-HEM, waste water (WasteWatR TM)**

Art. ID                      ERA-519QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 5 mL flame-sealed ampule yields up to 2 liters after dilution. Contains both hexadecane and stearic acid. Note: If a NELAC compliant PT is required, use Cat. #582 or Cat. #4120.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	HEM		5 - 100	mg/L		
	SGT-HEM		5 - 100	mg/L		