

## **QuiKResponse - Complex Nutrients, waste water (WasteWatR TM)**

Art. ID ERA-525QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Kjeldahl-Nitrogen as N		1,5 - 35	mg/L		
	Total Phosphorus as P		0,5 - 10	mg/L		