

QuiKResponse - Inorganic Disinfection 2

Art. ID ERA-5262QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 23 mL screw-cap vial yields up to 4 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bromate		7,0 - 50,0	µg/L		
	Bromide		50 - 300	µg/L		