

QuiKResponse - Inorganic Disinfection 1

Art. ID ERA-5272QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 23 mL screw-cap vial yields up to 4 liters after dilution. The diluted standard will contain chlorate at 60 to 180 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chlorate		60 - 180	µg/L		
	Chlorite		100 - 1000	µg/L		