

## **QuiKResponse - Oil & Grease (O&G) in Soil**

Art. ID	ERA-549QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One screw-cap bottle containing 50 g of soil ready to analyze. Use with gravimetric method 9071B or infrared spectrometric analysis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Hexane Extractable Material (O&G) (Gravimetric)		300 - 3000	mg/kg		
	n-Hexane Extractable Material (O&G) (IR)		300 - 3000	mg/kg		