

## **QuiKResponse - Oil & Grease (O&G) in Soil**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | ERA-549QR                         |
| Unit            | each                              |
| Deliverydetails | No Dangerous Good /not restricted |

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One screw-cap bottle containing 50 g of soil ready to analyze. Use with gravimetric method 9071B or infrared spectrometric analysis.

| Text/Information | Analyte/Parameter                                 | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|---|------------|---------------------|-------|--------|--------|
|                  | n-Hexane Extractable Material (O&G) (Gravimetric) |            | 300 - 3000          | mg/kg |        |        |
|                  | n-Hexane Extractable Material (O&G) (IR)          |            | 300 - 3000          | mg/kg |        |        |