

QuiKResponse - TPH in Soil, QuiK

Art. ID	ERA-572QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One screw-top bottle contains 50 g of soil with TPH. Use with EPA methods 8440 and 9071B. TPH is present in the range of 250-3,000 mg/kg.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	non-Polar Extractable Material (TPH) (Gravimetric)		300 - 3000	mg/kg		
	non-Polar Extractable Material (TPH) (IR)		300 - 3000	mg/kg		