

QuiKResponse - TPH in Soil, QuiK

| | |
|-----------------|-----------------------------------|
| Art. ID | ERA-572QR |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One screw-top bottle contains 50 g of soil with TPH. Use with EPA methods 8440 and 9071B. TPH is present in the range of 250-3,000 mg/kg.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|------------|---------------------|-------|--------|--------|
| | non-Polar Extractable Material (TPH) (Gravimetric) | | 300 - 3000 | mg/kg | | |
| | non-Polar Extractable Material (TPH) (IR) | | 300 - 3000 | mg/kg | | |