

QuiKResponse - TPH in Water, QuiK

Art. ID	ERA-602QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! This is sold as a QuiK? Response product only. One liter Whole-Volume bottle is ready to analyze for Total Petroleum Hydrocarbons. Use with EPA methods 413.1, 413.2, 418.1, 1664 and 5520.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Petroleum Hydroca rbons (Gravimetric)		20 - 170	mg/L		
	Total Petroleum Hydroca rbons (Infrared)		20 - 170	mg/L		