

QuiKResponse - Hexavalent Chromium, potable water (PotableWatR TM)

Art. ID ERA-658QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|--------------|---------------------|------|--------|--------|
| | Hexavalent chromium | [18540-29-9] | 5,0 - 50 | µg/L | | |