

## **QuiKResponse - Hexavalent Chromium, potable water (PotableWatR TM)**

Art. ID                      ERA-658QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hexavalent chromium	[18540-29-9]	5,0 - 50	µg/L		