

QuiKResponse - Vanadium, potable water (PotableWatR TM)

Art. ID ERA-660QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution. Designed to meet California ELAP requirements.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Vanadium (V)	[7440-62-2]	5,0 - 50	µg/L		