

QuiKResponse - UV 254 Absorbance, potable water (PotableWatR TM)

Art. ID ERA-662QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	UV 254 Absorbance		0,02 - 0,7	cm-1		