

## **QuiKResponse - o-Phosphate Nutrients, potable water (PotableWatR TM)**

Art. ID                      ERA-667QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

| Text/Information | Analyte/Parameter    | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------|------------|---------------------|------|--------|--------|
|                  | Ortho-Phosphate as P |            | 0,5 - 5,5           | mg/L |        |        |