

QuiKResponse - Chloral Hydrate, potable water (PotableWatR TM)

Art. ID ERA-676QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Includes chloral hydrate at 4-30 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chloralhydrate	[302-17-0]	4,0 - 30,0	µg/L		