

## **QuiKResponse - Haloacetic Acids (HAA), potable water (PotableWatR TM)**

Art. ID ERA-684QR  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2398 II E2

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Includes all the analytes below at 5-50 µg/L.

| Text/Information | Analyte/Parameter      | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|-------------|---------------------|------|--------|--------|
|                  | Bromochloroacetic acid | [5589-96-8] | 5,0 - 50,0          | µg/L |        |        |
|                  | Dibromoacetic acid     | [631-64-1]  | 5,0 - 50,0          | µg/L |        |        |
|                  | Dichloroacetic acid    | [79-43-6]   | 5,0 - 50,0          | µg/L |        |        |
|                  | Bromoacetic acid       | [79-08-3]   | 5,0 - 50,0          | µg/L |        |        |
|                  | Chloroacetic acid      | [79-11-8]   | 10,0 - 50,0         | µg/L |        |        |
|                  | Trichloroacetic acid   | [76-03-9]   | 5,0 - 50,0          | µg/L |        |        |