

QuiKResponse - Semivolatiles #2 Herbicides potable water (PotableWatR TM)

Art. ID ERA-691QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains all the analytes below at 8-800 µg/L.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Diquat | [85-00-7] | 8,0 - 40,0 | µg/L | | |
| | Endothal | [145-73-3] | 90 - 500 | µg/L | | |
| | Glyphosate | [1071-83-6] | 375 - 800 | µg/L | | |
| | Paraquat | [1910-42-5] | 8 - 100 | µg/L | | |