

## **QuiKResponse - Hardness, potable water (PotableWatR TM)**

Art. ID ERA-693QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 250 mL whole-volume bottle is ready to analyze.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	30 - 90	mg/L		
	Calcium hardness as CaC O3		75 - 225	mg/L		
	Magnesium (Mg)	[7439-95-4]	2,0 - 20,0	mg/L		
	Sodium (Na)	[7440-23-5]	12,0 - 50,0	mg/L		
	Total hardness as CaCO3		83 - 307	mg/L		