

### **QuiKResponse - Nitrite, potable water (PotableWatR TM)**

Art. ID                      ERA-695QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

#### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nitrite as N		0,4 - 2	mg/L		