

## **QuiKResponse - Residual Chlorine, potable water (PotableWatR TM)**

Art. ID                      ERA-696QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Residual free chlorine		0,5 - 3	mg/L		
	Total residual chlorine		0,5 - 3	mg/L		