

## **QuiKResponse - Inorganics, potable water (PotableWatR TM)**

Art. ID	ERA-698QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 500 mL whole-volume bottle is ready to analyze. The CRM is also certified for Sodium. For a Sodium PT, order Hardness, ERA-555.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alkalinity as CaCO <sub>3</sub>		25 - 200	mg/L		
	Chloride (Cl)		20 - 160	mg/L		
	Fluoride (F)		1 - 8	mg/L		
	Nitrate + Nitrite as N		3 - 10	mg/L		
	Nitrate (NO <sub>3</sub> -) as N		3 - 10	mg/L		
	Potassium (K)	[7440-09-7]	10 - 40	mg/L		
	Sodium (Na)	[7440-23-5]	10 - 400	mg/L		
	Conductivity (25°C)		130 - 1300	µmhos/cm		
	Sulfate		25 - 250	mg/L		
	Total Dissolved Solids at 180 °C (TDS)		100 - 1000	mg/L		