

## **QuiKResponse - Toxaphene, potable water (PotableWatR TM)**

Art. ID	ERA-700QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains toxaphene at 2-20 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Campechlor (Toxaphene) (total)		3,0 -20,0	µg/L		