

QuiKResponse - Halomethanes (THMs), potable water (PotableWatR TM)

Art. ID ERA-702QR
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains all of the analytes below at 5-50 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bromodichloromethane	[75-27-4]	5,0 - 50,0	µg/L		
	Bromoform (Tribromomethane)	[75-25-2]	5,0 - 50,0	µg/L		
	Chlorodibromomethane	[124-48-1]	5,0 - 50,0	µg/L		
	Chloroform (Trichloromethane)	[67-66-3]	5,0 - 50,0	µg/L		