

## **QuiKResponse - Chlorinated Acid Herbicides, potable water (PotableWatR TM)**

Art. ID                      ERA-704QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Includes at least 12 of the analytes below at 1-120 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,5-T	[93-76-5]	10 - 100	µg/L		
	2,4,5-TP	[93-72-1]	10 - 100	µg/L		
	2,4-D	[94-75-7]	10 - 100	µg/L		
	2,4-DB	[94-82-6]	20 - 120	µg/L		
	3,5-Dichlorobenzoic acid	[51-36-5]	10 - 100	µg/L		
	4-Nitrophenol	[100-02-7]	10 - 100	µg/L		
	Acifluorfen	[50594-66-6]	10 - 100	µg/L		
	Bentazone	[25057-89-0]	10 - 100	µg/L		
	Chloramben	[133-90-4]	20 - 100	µg/L		
	DCPA	[1861-32-1]	20 - 100	µg/L		
	Dalapon	[75-99-0]	10 - 100	µg/L		
	Dicamba	[1918-00-9]	20 - 100	µg/L		
	Dichlorprop	[120-36-5]	10 - 100	µg/L		
	Dinoseb	[88-85-7]	7,0 - 70	µg/L		
	Pentachlorophenol	[87-86-5]	1,0 - 25	µg/L		
	Picloram	[1918-02-1]	10 - 100	µg/L		