

QuiKResponse - Chlordane, potable water (PotableWatR TM)

Art. ID ERA-705QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains technical chlordane at 2-20 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chlordane (total)		2,0 - 20,0	µg/L		