

## **QuiKResponse - EDB/DBCP/TCP, potable water (PotableWatR TM)**

Art. ID ERA-706QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Each lot contains all analytes below at 0.05-2 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3-Trichloropropane	[96-18-4]	0,1 - 2	µg/L		
	1,2-Dibromo-3-chloropropane	[96-12-8]	0,1 - 2	µg/L		
	Ethylene Dibromide	[106-93-4]	0,1 - 2	µg/L		