

## **QuiKResponse - PCBs as Decachlorobiphenyl, potable water (PotableWatR TM)**

Art. ID	ERA-708QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Includes an Aroclor randomly selected from the list below at 0.5-5 µg/L as decachlorobiphenyl.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aroclor 1016	[12674-11-2]	0,26 - 2,6	µg/L		
	Aroclor 1221	[11104-28-2]	0,19 - 1,9	µg/L		
	Aroclor 1232	[11141-16-5]	0,23 - 2,3	µg/L		
	Aroclor 1242	[53469-21-9]	0,26 - 2,6	µg/L		
	Aroclor 1248	[12672-29-6]	0,3 - 3,0	µg/L		
	Aroclor 1254	[11097-69-1]	0,33 - 3,3	µg/L		
	Aroclor 1260	[11096-82-5]	0,36 - 3,6	µg/L		
	PCBs as decachlorobiphenyl		0,5 - 5	µg/L		