

QuiKResponse - Pesticides, potable water (PotableWatR TM)

Art. ID	ERA-709QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains at least 16 analytes randomly selected from the list below at 0.2-20 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alachlor	[15972-60-8]	2,0 - 20,0	µg/L		
	Aldrin	[309-00-2]	0,4 - 2	µg/L		
	Atrazine	[1912-24-9]	3,0 - 30,0	µg/L		
	Bromacil	[314-40-9]	2,0 - 20,0	µg/L		
	Butachlor	[23184-66-9]	8,0 - 80,0	µg/L		
	Diazinon	[333-41-5]	0,1 - 100	µg/L		
	Dieldrin	[60-57-1]	0,5 - 3	µg/L		
	Endrin	[72-20-8]	0,1 - 5	µg/L		
	Heptachlor	[76-44-8]	0,4 - 5	µg/L		
	Heptachlor Epoxide (beta)	[1024-57-3]	0,2 - 5	µg/L		
	Hexachlorobenzene	[118-74-1]	0,5 - 4	µg/L		
	Hexachlorocyclopentadiene	[77-47-4]	2,0 - 30,0	µg/L		
	Methoxychlor	[72-43-5]	10 - 100	µg/L		
	Metolachlor	[51218-45-2]	8,0 - 80,0	µg/L		
	Metribuzin	[21087-64-9]	2,0 - 60,0	µg/L		
	Molinate (Ordram)	[2212-67-1]	5,0 - 50	µg/L		
	Prometon	[1610-18-0]	0,1 - 100	µg/L		
	Propachlor	[1918-16-7]	1,0 - 4,0	µg/L		
	Simazine	[122-34-9]	4,0 - 40,0	µg/L		
	Thiobencarb	[28249-77-6]	0,1 - 100	µg/L		
	Trifluralin	[1582-09-8]	1,0 - 5,0	µg/L		
	gamma-HCH (Lindane)	[58-89-9]	0,2 - 5	µg/L		