

## **QuiKResponse - Acids, priority pollutant (PriorityPollutnT TM)**

Art. ID ERA-712QR  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Contains a subset of the analytes listed below at 30-200 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,4,6-Tetrachlorophenol	[58-90-2]	30 - 200	µg/L		
	2,4,5-Trichlorophenol	[95-95-4]	30 - 200	µg/L		
	2,4,6-Trichlorophenol	[88-06-2]	30 - 200	µg/L		
	2,4-Dichlorophenol	[120-83-2]	30 - 200	µg/L		
	2,4-Dimethylphenol	[105-67-9]	40 - 200	µg/L		
	2,4-Dinitrophenol	[51-28-5]	100 - 200	µg/L		
	2,6-Dichlorophenol	[87-65-0]	30 - 200	µg/L		
	2-Chlorophenol	[95-57-8]	30 - 200	µg/L		
	2-Methylphenol	[95-48-7]	40 - 200	µg/L		
	2-Nitrophenol	[88-75-5]	50 - 200	µg/L		
	4,6-Dinitro-2-methylphenol (DNOC)	[534-52-1]	40 - 200	µg/L		
	4-Chloro-3-methylphenol	[59-50-7]	30 - 200	µg/L		
	4-Methylphenol	[106-44-5]	50 - 200	µg/L		
	4-Nitrophenol	[100-02-7]	100 - 200	µg/L		
	Benzoic acid	[65-85-0]	30 - 200	µg/L		
	Pentachlorophenol	[87-86-5]	40 - 200	µg/L		
	Phenol	[108-95-2]	100 - 200	µg/L		