

## **QuiKResponse - Chlordane, priority pollutant (PriorityPollutnT TM)**

Art. ID ERA-716QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters of sample after dilution. Contains technical chlordane at 3-25 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chlordane (total)		3,0 - 25,0	µg/L		