

## **QuiKResponse - Toxaphene, priority pollutant (PriorityPollutnT TM)**

Art. ID ERA-717QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters of sample after dilution. Contains toxaphene at 20-100 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Campechlor (Toxaphene)	[8001-35-2]	20 - 100	µg/L		