

QuiKResponse - BTEX & MTBE in Water

Art. ID	ERA-760QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Includes all BT EX compounds and MTBE at 7-300 µg/L after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	8 - 120	µg/L		
	Ethylbenzene	[100-41-4]	15 - 100	µg/L		
	Methyl-tert.butylether (MTBE)	[1634-04-4]	9 - 100	µg/L		
	Toluene	[108-88-3]	7 - 100	µg/L		
	Xylenes, total	[108-38-3]/[106 -42-3]/[95-47-6]	20 - 300	µg/L		