

QuiKResponse - Gasoline Range Organics (GRO) in Water

| | |
|-----------------|-----------------------------------|
| Art. ID | ERA-762QR |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with both purge & trap and modified EPA 8015 GC/FID methods to test for GRO at 400-4,000 µg/L. Also use to test for BT EX in gasoline.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------|------------|---------------------|------|--------|--------|
| | Benzene in GRO | [71-43-2] | 0,5 - 1000 | µg/L | | |
| | Ethylbenzene in GRO | [100-41-4] | 0,5 - 1000 | µg/L | | |
| | Gasoline Range Organics | | 200 - 4000 | µg/L | | |
| | Toluene in GRO | [108-88-3] | 0,5 - 1000 | µg/L | | |
| | Xylenes, total in GRO | | 20 - 450 | µg/L | | |