

QuiKResponse - Gasoline Range Organics (GRO) in Water

Art. ID	ERA-762QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with both purge & trap and modified EPA 8015 GC/FID methods to test for GRO at 400-4,000 µg/L. Also use to test for BT EX in gasoline.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene in GRO	[71-43-2]	0,5 - 1000	µg/L		
	Ethylbenzene in GRO	[100-41-4]	0,5 - 1000	µg/L		
	Gasoline Range Organics		200 - 4000	µg/L		
	Toluene in GRO	[108-88-3]	0,5 - 1000	µg/L		
	Xylenes, total in GRO		20 - 450	µg/L		