

## **QuiKResponse - Diesel Range Organics (DRO) in Water**

Art. ID ERA-764QR  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Includes #2 Diesel at 800-6,000 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diesel Range Organics		800 - 6000	µg/L		