

QuiKResponse - Surfactants-MBAS, waste water (WasteWatR TM)

Art. ID ERA-776QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Surfactants (MBAS)		0,2 - 1	mg/L		