

QuiKResponse - Surfactants-MBAS, waste water (WasteWatR TM)

Art. ID ERA-776QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|------|--------|--------|
| | Surfactants (MBAS) | | 0,2 - 1 | mg/L | | |