

QuiKResponse - Turbidity, waste water (WasteWatR TM)

Art. ID ERA-777QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 24 mL screw-cap vial yields up to 1 liter after dilution. Use with nephelometric methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Turbidity		2,0 - 30,0	NTU		