

## **QuiKResponse - Turbidity, waste water (WasteWatR TM)**

Art. ID                      ERA-777QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 24 mL screw-cap vial yields up to 1 liter after dilution. Use with nephelometric methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Turbidity		2,0 - 30,0	NTU		