

## **QuiKResponse - Surfactants-MBAS, potable water (PotableWatR TM)**

Art. ID                      ERA-784QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	MBAS		0,05 - 1	mg/L		