

## **QuiKResponse - Texas Low-Level Fuels (TPH) in Water**

Art. ID                      ERA-794QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains unleaded regular gasoline and No. 2 Diesel Fuel resulting in TPH in the range 5-10 mg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total petroleum hydrocarbons (TPH)		5,0 - 10,0	mg/L		