

QuiKResponse - Gasoline Additives, potable water (PotableWatR TM)

Art. ID ERA-909QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains all of the analytes below at 5-50 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diisopropylether	[108-20-3]	5,0 - 50	µg/L		
	Ethyl-tert-butylether (ETBE)	[637-92-3]	5,0 - 50	µg/L		
	Methyl-tert.butylether (MTBE)	[1634-04-4]	5,0 - 50	µg/L		
	tert.-Amylmethylether (TAME)	[994-05-8]	5,0 - 50	µg/L		
	Trichlorotrifluoroethane		5,0 - 50	µg/L		
	tert.-Butanol	[75-65-0]	5,0 - 50	µg/L		