

## **QuiKResponse - Perchlorate, potable water (PotableWatR TM)**

Art. ID ERA-910QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Perchlorate		4,0 - 20,0	µg/L		