

QuiKResponse - Settleable Solids, waste water (WasteWatR TM)

Art. ID ERA-911QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 60 mL poly bottle with a solid yields 1 liter after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Settleable Solids | | 5 - 100 | mL/L | | |