

### **QuiKResponse - Volatile Solids, waste water (WasteWatR TM)**

Art. ID                      ERA-913QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

#### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 12 mL screw-cap vial with a solid yields 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Volatile Solids		100 - 500	mg/L		