

## **QuiKResponse - Acidity, waste water (WasteWatR TM)**

Art. ID ERA-915QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 250 mL whole-volume bottle is ready to analyze. Designed for use with titrimetric methods to a pH endpoint of 8.3.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acidity as CaCO <sub>3</sub>		650 - 1800	mg/L		