

QuiKResponse - Acidity, waste water (WasteWatR TM)

Art. ID ERA-915QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 250 mL whole-volume bottle is ready to analyze. Designed for use with titrimetric methods to a pH endpoint of 8.3.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acidity as CaCO ₃		650 - 1800	mg/L		