

QuiKResponse - Low-Level Total Residual Chlorine

Art. ID ERA-917QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! Designed for testing at low µg/L levels. One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with sensitive titrimetric or colorimetric methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total residual chlorine		20 - 250	µg/L		