

QuiKResponse - Uranium, potable water (PotableWatR TM)

Art. ID ERA-930QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 500 mL whole-volume bottle is ready to analyze for corrosivity, calcium carbonate saturation and Langelier saturation index.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Uranium (U)	[7440-61-1]	3 - 104	µg/L		